

NB200

Strip Pouch Packager

NORITSU

Pharmacy Automation

Automatically Fills and Validates Multidose Blister Cards


200 hot-swappable
cassettes

In-line automated
validation

Process 35 cards/hour *

Handles half-tabs and
infrequent medications

Multi-unit scalability



Powered By



NexusRx is a **comprehensive workflow management system** that integrates with all Noritsu automation equipment, managing everything from prescription preparation and production to reporting and verification.

Get a complete, streamlined view of your operations, while optimizing your workflow. Coordinate dispensing, tracking medications, and reporting from anywhere in the facility. Achieve higher accuracy rates, better manage labor, and reduce operational bottlenecks with NexusRx.


Streamline your blister card checking and verification process with the **BCO** workstation.

- Visually check the contents of any card and quickly scan for any flagged issues
- Add controlled dispenses, such as Cytotoxic or hormone therapy medications manually
- Manually fill and check one-off blister cards



Specifications

Example floor lay-out for blister packing with the



Measurements	76.8" x 33.5" x 79.1" (W x D x H)
Required footprint	8.2' x 4.9' (W x D)
Weight	1,808 lbs (820 kg)
Electrical parameters	110 ~ 115 Vac - 50/60 Hz - 7.3A (800VA)
Network requirements	1 Gbit/s minimum network connection
Operating environment	41 ~ 104 °F and < 80 %RH (non-condensing)
Maximum card size	10.4" x 8.7" x 1"
Pill size limitations	Length/width 5 - 23 mm Thickness 1.5 - 8.5 mm
Maximum throughput	Up to 40 trays/hour
Noise level	58 ~ 63 dB(A)

Noritsu/GF-VBM-2302